

VWR International Ltd
1-3 The Markham Centre
Station Road
Theale
Reading
RG7 4PE

23rd October 2020

Discontinuation: Avantor® ACE® Excel 3 µm & 5 µm columns for CN & Sil chemistries

Dear Chromatographer

In an effort to streamline our product offering, we will no longer be offering two options for the Avantor® ACE® CN and Sil chemistries with 3 µm and 5 µm particle size. Effective immediately, we will offer only the ACE- HPLC pressure rated (275 bar) options for these two stationary phases (with the exception of 250mm length, 3 µm which generate high back pressures and remain in EXL- UHPLC compatible hardware). The EXL- high pressure (1,000 bar) versions will be discontinued after any remaining stocks have been exhausted.

Included with this letter is a list of those catalogue numbers affected by this change alongside the recommended replacement column which contains identical packing material (particle size, chemistry etc.) but is pressure rated to HPLC pressures.

We value any feedback you have, please let me know if you have any questions or feel free to contact our Technical support team : chromsupport@avantorsciences.com

Kind Regards



Gemma Lo

Senior Product Specialist - Chromatography